

Kinetec Fixed Treatment Table



Blocking bar on the foot end of the table



Head rest (optional)



Showned with Kinetec Spectra Knee CPM Machine

Our treatment table is **especially** adapted to be used in combined in conjunction with our rehabilitation devices.



The blocking bar prevents damage to the upholstery and maintains a comfortable position for the patient.

Large width of the tabletop (74cm – 29 Inch) improves the stability of the machine and comfort for the patient, there is sufficient space for the non-mobilised leg on the table.



Lowest table (55cm – 22 Inch) facilitating access to patients in wheelchairs or with reduced mobility.

Tabletop of 2,05m (81 Inch) with backrest of 72cm (28 Inch) allows tall patients to be comfortably installed, the backrest is adjustable from 0 ° to 60 °. For the oversized patients a head rest is available as an optional accessory.

Comfortable upholstery, foam thickness of 60mm (2.4 Inch) high impact strength density of 40Kg (88lb) with treated cover PVC M1.

Steel fixed frame allows patients up to 180Kg (397lb).

The Kinetec Fixed Treatment Table is available in the colours here below shown, (need to be stipulated when placing an order), 5 weeks delivery time for all colours except for the colour Grey* in standard configuration.



Grey* Blue Dark bue Dark grey Light green Dark Green Red

Technical Specifications:

- Maximum user weight: 180Kg (397lb)
- Upholstery in PVC:
 - Thickness 1,1mm (0,43Inch)
 - Composition 84% PVC, 16% Coton
 - Fire Classification M1
- Total weight of the table: 24Kg
- Require assembling
(delivered in 3 parts: the tabletop, the base, the blocking bar)

Ordering Information

- 467500027 Kinetec Fixed Treatment Table
- 467500035 Head rest (optional)
- Complete rehabilitation station**
- 4621708042 Table with Prima Advance
- 4621706042 Table with Spectra
- 4621701842 Table with Performa



Kinetec Electric Treatment Table



Showned with Kinetec Prima Advance Knee CPM Machine

Our treatment table is **especially** adapted to the combined use with our rehabilitation devices.

The blocking bar prevents damage to the upholstery and maintains a comfortable position for the patient.

The blocking bar is adjustable in length allowing tall patients to be correctly installed.

Technical characteristics are identical to the 2 section electric table on page 23, available in the same colours.




Blocking bar on the table foot end

Ordering Information

4675000043 Kinetec Electric Treatment Table



The 2 and 3 section electric treatment tables are available in colours shown here below, which needs to be stipulated when placing an order, 5 weeks delivery time for all colours except for the colour Grey* in standard configuration.




With Grey* Beige Black Dark Blue Light Green Lilac Blue Pink Red Yellow Sky Blue

Technical Specifications:

- Manual hand set
- Current Consumption: 0.9A
- Power: 175W
- Protective Class 2.B
- Degree of Protection IP66
- Mains Plug Fuse Rating: 5A

- Vinyls comply with BS5852, Part1; Ignition Sources 0&1 (When tested in combination with polyurethane foam) and BS5852, Part2; Ignition Source 5 (when tested in combination with CMHR foams)
- Warranty: 7 years on the frame
 3 years on the motor
 1 year on the upholstery

Electric Treatment Table



Can be used in hospitals, medical and rehabilitation centres as well as in health centres and private practices for massage and manual therapy.

- **Maximum user weight: 180Kg (397lb)**
- **Luxury Upholstery**
- **Variable height from 49 to 94 cm (19 to 37 Inch)**
- **Standard retractable castors of 50mm (2 Inch)**
- **Adjustable feet for uneven floors**
- **Wide range of accessories**



Ordering Information

| | |
|------------|------------------------------------|
| 4665025002 | 2 section Electric Treatment Table |
| 4665035002 | 3 section Electric Treatment Table |

The treatment table provides great value with strength, elegance and reliability and includes:

- Breathing Hole, Chrome Hanging Rail, Gas Springs to assist the backrest
- Grey: the standard colour for fast delivery, (10 other colours available, see page 22)

Accessories for 2 and 3 section Electric Treatment Table:

- 1 4665004790 Paper Roll Holder
- 2 4665004873 Adjustable Arm Rest
- 3 4665004881 Arm Support Tables (Upholstered)
- 4 4665004823 Painted Cot Sides (not illustrated)
- 5 4665004831 Painted Cot Sides Padded (not illustrated)
- 6 4665004849 Chromed Cot Sides
- 7 4665004857 Chromed Cot Sides Padded
- 8 4665004865 Cordless Remote Control
- 9 4665015402 Foot Switch
- 10 4665015202 Electric Backrest
- 11 4665015302 Electric Tilting Seat (not illustrated) (for 3 sections table only)
- 12 4665004815 Transportation Castors \varnothing 100mm (4Inch)

